

LF8.001.000 [LabFlame basic](#), complete with full stainless steel footswitch



= Basic unit including stainless steel footswitch and all LabFlame advantages such as stainless steel housing and:

- 3 standard-programs for foot pedal and button (function knob)
- SCS (Safety Control System), BHC (Burner Head Control)

LF8.003.000 [LabFlame RF](#), with wireless radio footswitch, stainless steel



- Same range as LabFlame basic
- with integrated radio system and wireless foot pedal (stainless steel)

LF8.000.000 [LabFlame IR](#), safety burner with IR-Sensor and LC Display



= Enhanced safety burner with IR sensor
Offers all LabFlame basic advantages and:

- IR-Sensor and LC-Display
- Conventional and exclusive DoubleClick IR-Sensor
- Adjustable IR-Sensor reaction distance
- 5 standard-programs for IR-Sensor, button (function knob)
- Holding device for 3 inoculation loop holders
- Optional foot pedal see accessory price list

LF8.004.000 [LabFlame +](#) with IR-Sensor, LCD + integrated NIMI (AA) Rechargeable battery pack build in



= Enhanced safety burner with IR sensor
Offers all LabFlame IR advantages and cable free operation
Based on NIMI rechargeable battery pack build in



[Ready to go](#) all LabFlame units do come fully equipped and are ready to go

- Removable and decomposable burner head optional life science burner head
- Tilt mechanism right / left
- Holding device for 3 inoculation loop holders
- Nozzles for natural, propane/butane gas
- Turbo flame
- Wrench 17mm for gas connection
- Screwdriver for burner head and cover of burner shaft
- Tubing connector with swivel nut, switching power supply (global)